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The ROBOT DOSER is a Tube-free design Laboratory dosing system for the ultimate precision in dye bath preparation and allow textile laboratory to improve reproducibility and efficiency without adding personnel.

The functioning principle is based on gravimetric measuring by means of robot pipette and electronic balance. After input of the recipe, robot pipette moves to the dyestuff solution and sucks the required dyestuff liquid, then inject to the pots of dyeing machine and weighed gravimetrically by electronic balance.

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Using high technology robot pipette and electronic balance weighing system, to ensure the accuracy of liquid dispensing, no more problem of unwanted drip contamination that often occurs in tube design system, now all tubes, solenoid valve troubles, such as sedimentation of dyes, are history in our ROBOT DOSER.

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
Thanks to the independent injector design, each bottles equipped with own injector, no accidental contamination risks of the various solutions in between, furthermore, traditional intermediate washing, drying cycles of injector is not required, thus increasing the working efficiency in comparison with single injector system, no any lost of time for cleaning the injector.

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
Stainless steel made injector to replace traditional glass or plastic material injector.

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
Gravimetric measuring system by means of Swiss made electronic scale, not traditional volumetric machine, the weight of each solution dispensed can be controlled in real time, the degree of accuracy achieved is without any doubt. An external monitor connected to electronic scale for easy observation of on-line weighing status.

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
Modular structure made of aluminium and coated steel, windshield design a protection panel completely isolates the working area from the outside, ensuring perfect environmental conditions and safety.

可连续式滴液 / Continuous weighing capability :

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ROBOT DOSER offers a higher level of automation, the robotized arm inject the solutions into dye pot where 1 unique conveyor system holding up to 30 dye pots with continuous weighing capability.

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- ◆ Models available in 120 bottle is (or 80 bottles), each bottle is equipped with independent injector
- ◆ Robotic X-Y-Z axis mechanism, complete with 3 servo motors
- ◆ Precise electronic balance with 0.01 g readability
- ◆ 5 dyepot trays for dispensing recipes (with unique conveyor system holding up to 30 dyepots)
- ◆ Dispense accuracy of 0.02 grams
- ◆ Automatic solution bottle agitation, with speed control
- ◆ Complete suite of Dispensing and management software
- ◆ Dimension (L x W x H) :
319 x 146 x 198 cm (Model RD-120)
260 x 146 x 198 cm (Model RD-80)

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/Unique conveyor system

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/Gravimetric dispensing

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/Dispensing of water up to 6 pot simultaneously

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/Aluminum frame + cover the steel board

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Laborteks Analitik ve Kalite Kontrol Test Cihazları San.Tic.Ltd.şti